

Cross straps



FIELD OF APPLICATION

Removable and adjustable cross straps for metallic aprons.

MATERIAL : Polyurethane

FEATURES

Colour :Blue

Thickness : 1000 µm



Ref. OATP010.50.000

CLEANING

The cross straps should undergo the same cleaning process as the metallic apron. You should rinse the cross straps with a minimum 80°C water temperature. You should disinfect the cross straps with a non-acid product authorized by the health services.

Do not use any chlorine or chlorinated alkaline products at high concentration, and avoid prolonged contact with these products

N° intrastat : 39262000

Bien plus que la sécurité
More than just safety