

## VITOFIX harness



### FIELD OF APPLICATION

Removable and adjustable harness for metallic aprons.

**MATERIAL:** Polyurethane


### FEATURES

**Colour:** blue

**Thickness:** 1 000 µm

### SIZE AND REFERENCES

Reference	Size	Chest measurement (cm)
0VVITO.001.05	Size 1	< 100 cm
0VVITO.002.05	Size 2	from 90 to 120 cm
0VVITO.003.05	Size 3	> 120 cm



Size 1

Size 2

Size 3

### CLEANING

The VITOFIX harness should undergo the same cleaning process as the metallic apron. You should disinfect the VITOFIX harness with a non-acid product authorized by the health services. You should rinse the VITOFIX harness with a minimum 80°C water temperature.

Do not use any chlorine or chlorinated alkaline products at high concentration, and avoid prolonged contact with these products.

N° intrastat : 39262000