


VITOFIX HARNESS



MATERIAL

Polyurethane

FEATURES

- Colour : blue
- Thickness : 1000 µm

FIELD OF APPLICATION

Removable and adjustable harness for metallic aprons.


CLEANING CONDITIONS

The VITOFIX harness should undergo the same cleaning process as the metallic apron. You should disinfect the VITOFIX harness with a non-acid product authorized by the health services. You should rinse the VITOFIX harness with a minimum of 80°C water temperature. Do not use any chlorine or chlorinated alkaline products at high concentration, and avoid prolonged contact with these products.

SIZES AND REFERENCES

Reference	Size	Chest measurement (in cm)
0VVITO.001.05	Size 1	< 100 cm
0VVITO.002.05	Size 2	from 90 to 120 cm
0VVITO.003.05	Size 3	> à 120 cm

SIZE 1



SIZE 2



SIZE 3

